Abbott

NycoCard™
CRP

EN	5-8	FI	29-32	NO	53-56
CS	9-12	FR	33-36	PL	57-60
DA	13-16	HU	37-40	PT	61-64
DE	17-20	IT	41-44	RO	65-68
EL	21-24	LT	45-48	RU	69-72
ES	25-28	NL	49-52	SV	73-76

REF 1116807

1116924 Rev. A 2019/02

SYMBOLS/ SYMBOLY/ SYMBOLER/ SYMBOLE/ ΣΥΜΒΟΛΑ/ SÍMBOLOS/ SYMBOLIT/ SYMBOLES/ JELMAGYARÁZAT/ SIMBOLI/ SIMBOLIAI/ SYMBOLEN/ SYMBOLER/ SYMBOLE/ SÍMBOLO/ SIMBOLURI/ СИМВОЛЫ/ SYMBOLER

Conformity to the European directive 98/79/EC on *in vitro* diagnostic medical devices
 Splnění evropské směrnice 98/79/ES o zdravotnických prostředcích pro diagnostiku *in vitro*
 Overensstemmelse med EU-direktivet 98/79/EC om *in vitro* diagnostisk udstyr til medicinsk brug
 Konformität mit der europäischen Richtlinie 98/79/EC über medizinische Geräte zur *In-vitro*-Diagnostik
 Συμμορφώνεται με την ευρωπαϊκή οδηγία 98/79/ΕΚ περί ιατροτεχνολογικών βοηθημάτων που χρησιμοποιούνται στη διάγνωση *in vitro*
 Conformidad con la Directiva europea 98/79/CE sobre dispositivos médicos para el diagnóstico *in vitro*
 Yhdenmukainen Euroopan direktiivin 98/79/EC kanssa koskien *in vitro* diagnostisia lääkinnällisiä laitteita
 Conformité avec la directive européenne 98/79/CE sur les appareils médicaux de diagnostic *in vitro*
 Az Európai Parlament és a Tanács 98/79/EK irányelve az *in vitro* diagnosztikai orvostechnikai eszközökről
 Conformità con la direttiva europea 98/79/CE sui dispositivi medicali per la diagnostica *in vitro*
 Atitiktis Europos direktyvai 98/79/EB
 Voldoet aan de Europese richtlijn 98/79/EG betreffende medische hulpmiddelen voor *in-vitro*diagnostiek
 Overensstemmelse med EU-direktiv 98/79/EC for *in vitro* diagnostisk medisinsk utstyr
 Zgodność z dyrektywą 98/79/WE w sprawie wyrobów medycznych używanych do diagnostyki *in vitro*
 Conformidade com a Diretiva Europeia 98/79/CE de dispositivos médicos de diagnóstico *in vitro*
 Conformitate cu directiva europeană 98/79/CE privind dispozitivele medicale pentru diagnostic *in vitro*
 Знак соответствия европейской директиве 98/79/EC по медицинским устройствам, предназначенным для диагностики *in vitro*
 Överensstämmer EU-direktiv 98/79/EG om medicintekniska produkter för *in vitro*-diagnostik

Catalogue number
 Katalogové číslo
 Katalognummer
 Bestellnummer
 Αριθμός καταλόγου
 Número de catálogo
 Luettelonumero
 Numéro de référence du catalogue
 Katalógusszám
 Numero di catalogo
 Katalogo numeris
 Catalogusnummer
 Katalognummer
 Numer katalogowy
 Número do catálogo
 Număr de catalog
 Номер по каталогу
 Artikelnummer

Caution
 Varování
 Forsiktig
 Achtung
 Προσοχή
 Advertencia
 Varoitus
 Faire attention
 Figyelmeztetés
 Attenzione
 Įspėjimas
 Opgeet
 Advarsel
 Ostrzeżenie
 Atenção
 Avertizare
 Предупреждение
 Varning

In vitro diagnostic medical device
 Diagnostické lékařské zařízení pro testy *in vitro*
In vitro diagnostisk medicinsk utstyr
 Medizinisches Gerät zur *In-vitro*-Diagnostik
 Ιατροτεχνολογικό βοήθημα που χρησιμοποιείται στη διάγνωση *In vitro*
 Dispositivo médico para el diagnóstico *in vitro*
In vitro diagnostinen lääkinnällinen laite
 Dispositif médical de diagnostic *in vitro*
In vitro diagnosztikai orvosi eszköz
 Dispositivo medico per la diagnostica *in vitro*
In vitro diagnostinis medicinos prietaisai
 In-vitro-diagnostisch medisch apparaat
 Medisinsk utstyr til *in vitro* diagnostikk
 Wyrób medyczny do diagnostyki *in vitro*
 Dispositivo médico de diagnóstico *in vitro*
 Dispozitiv medical pentru diagnostic *in vitro*
 Медицинское устройство для диагностики *in vitro*
 Medicinteknisk produkt för *in vitro*-diagnostik

Lot number
 Číslo šarže
 Lot-nummer
 Chargen-Nr.
 Αριθμός παρτίδας
 Número de lote
 Eränumero
 Numéro de lot
 Gyártási tételszám
 Numero di lotto
 Partijos numeris
 Partijnummer
 Lotnummer
 Numer partii
 Número de lote
 Număr de lot
 Номер партии
 Lotnummer

2°C 8°C

Storage temperature 2-8°C (36-46°F)
Skladovací teplota 2-8 °C (36-46 °F)
Opbevaringstemperatur 2-8°C (36-46°F)
Lagertemperatur 2-8°C (36-46°F)
Θερμοκρασία φύλαξης 2-8°C (36-46°F)
Temperatura de almacenamiento 2-8°C (36-46°F)
Säilytyslämpötila 2-8°C (36-46°F)
Température de conservation 2-8°C (36-46°F)
Tárolási hőmérséklet 2-8 °C (36-46 °F)
Temperatura di conservazione 2-8°C (36-46°F)
Laikymo temperatūra (viršutinė ir apatinė riba)
Bewaartemperatuur 2-8°C (36-46°F)
Lagringstemperatur 2-8°C (36-46°F)
Temperatura przechowywania 2-8°C (36-46°F)
Temperatura de conservaço 2-8°C (36-46°F)
Temperatură de depozitare 2-8°C (36-46°F)
Температура хранения 2-8 °C (36-46 °F)
Förvaringstemperatur 2-8°C (36-46°F)

Consult the instruction for use
Viz návod k použití
Se brugervejledningen
Bitte die Gebrauchsanleitung beachten
Συμβουλευτείτε τις οδηγίες χρήσης
Consulte las instrucciones de uso
Katso lisätietoja käyttöohjeesta
Consulter le mode d'emploi
További információt a használati útmutatóban talál
Consultare le istruzioni per l'uso
Skaitykite naudojimo instrukciją
Raadpleeg de gebruiksaanwijzing
Se bruksanvisningen
Sprawdź w instrukcji użytkowania
Consulte as instruções de utilização
Consultati instrucțiunile de utilizare
См. инструкцию по использованию
Se bruksanvisningen

Contents sufficient for 48 tests
Obsah postačuje pro 48 testů.
Indholdet er tilstrækkeligt til 48 test
Inhalt ausreichend für 48 Tests
Περιεχόμενα επαρκή για 48 δοκιμές
Contenido suficiente para 48 pruebas
Sisältö riittää 48 testimäärälle
Contenu suffisant pour "48" dosages
48 vizsgálatra elegendő tartalom
Contenuto sufficiente per 48 test
Turinio pakanka 48 testų
Inhoud voldoende voor 48 testen
Innholdet rekker til 48 tester
"Zawiera ilość wystarczającą na 48 badań
Conteúdo suficiente para 48 testes
Conține o cantitate suficientă pentru 48 teste
Содержимого достаточно для проведения 48 анализов
Innehållet räcker till 48 analyser

Expiry date (year-month-day)
Datum expirace (rok-měsíc-den)
Udløbsdato (år-måned-dag)
Verfallsdatum (Jahr-Monat-Tag)
Ημερομηνία λήξης (έτος-μήνας-ημέρα)
Fecha de caducidad (año-mes-día)
Viimeinen käyttöpäivämäärä
(vuosi-kukausi-päivä)
Date de péremption (année-mois-jour)
Lejárati dátum (év-hónap-nap)
Data di scadenza (anno-mese-giorno)
Tinkamumo laikas (metai-mėnuo-diena)
Uiterste gebruiksdatum (jaar-maand-dag)
Utlepsdato (år-måned-dag)
Data przydatności (rok-miesiąc-dzień)
Prazo de validade (ano-mês-dia)
Data expirării (an-lună-ziuă)
Срок годности (год-месяц-день)
Utgångsdatum (år-månad-dag)

Date of Manufacture (year-month-day)
Datum výroby (rok-měsíc-den)
Produktionsdato (år-måned-dag)
Herstellungsdatum (Jahr-Monat-Tag)
Ημερομηνία κατασκευής
(έτος-μήνας-ημέρα)
Fecha de fabricación (año-mes-día)
Valmistuspäivämäärä
(vuosi-kukausi-päivä)
Date de fabrication (année-mois-jour)
Gyártás dátuma (év-hónap-nap)
Data di fabbricazione
(anno-mese-giorno)
Pagaminimo data (metai-mėnuo-diena)
Productiedatum (jaar-maand-dag)
Produksjonsdato (år-måned-dag)
Data produkcji (rok-miesiąc-dzień)
Data de fabrico (ano-mês-dia)
Data fabricației (an-lună-ziuă)
Дата изготовления (год-месяц-день)
Tillverkningsdatum (år-månad-dag)

Manufacturer
Výrobce
Producent
Hersteller
Κατασκευαστής
Fabricante
Valmistaja
Fabricant
Gyártó
Produttore
Gamintojas
Fabrikant
Produsent
Producent
Fabricante
Producător
Производитель
Tillverkare

For professional near-patient testing and laboratory use. For use with the NycoCard™ Reader II.

⚠ The NycoCard Reader II is calibrated for use with the NycoCard Tests only.

PRODUCT DESCRIPTION

Intended use

NycoCard™ CRP is an *in vitro* diagnostic test for the quantitative determination of C-reactive protein (CRP) in human whole blood and in human serum and plasma. The measurement of CRP provides information for the detection and evaluation of infection, tissue injury, inflammatory disorders and associated diseases.

Summary and explanation of the test

C-reactive protein (CRP) is one of the cytokine induced acute-phase proteins, the levels of which rise during a general, unspecific response to infections and non-infectious inflammatory processes^{3,5}. In healthy persons the serum or plasma CRP levels are below 5 mg/L¹. This threshold is often exceeded within four to eight hours after an acute inflammatory event, with CRP values reaching approximately 20 to 500 mg/L². As elevated CRP levels are always associated with pathological changes, the CRP assay provides information for the diagnosis, therapy, and monitoring of inflammatory diseases²⁻⁵.

Test principle

NycoCard CRP is a solid phase, sandwich-format, immunometric assay. In the test well of the device there is a membrane coated with immobilised CRP-specific monoclonal antibodies. A diluted sample is applied to the test device. When the sample flows through the membrane, the C-reactive proteins are captured by the antibodies. CRP trapped on the membrane will then bind the gold-antibody conjugate added, in a sandwich-type reaction. Unbound conjugate is removed from the membrane by the washing solution. A paper layer underneath the membrane absorbs excess liquid. In the presence of a pathological level of CRP in the sample, the membrane appears red-brown with colour intensity proportional to the CRP concentration of the sample. The colour intensity is measured quantitatively with the NycoCard Reader II.

Kit contents, 48 Tests

TD/Test Device	1 x 48 units
Plastic device containing a membrane coated with monoclonal anti-CRP antibodies.	
R1/Dilution Liquid	1 x 48 x 0.4 mL
Borate buffer (pH 9.0) and detergents.	
R2/Conjugate	1 x 3.5 mL
Solution containing monoclonal anti-CRP antibodies labeled with ultra-small gold particles.	
R3/Washing Solution	1 x 3.0 mL
Phosphate buffered NaCl solution (pH 7.4) and detergents.	

Material required

(not supplied with the kit)

- Pipette (50 µL) and pipette tips
- Capillary tube holder
- NycoCard Reader II
- Capillary tubes 5 µL
- Afinion™ CRP Control

Warnings and precautions

- For *in vitro* diagnostic use.
- Do not mix components from different kit containers.
- Do not use test kits after the expiry date.
- The reagents contain sodium azide in concentrations < 0.1% as a preservative. Sodium azide is a toxic agent. Avoid contact with eyes and skin. Flush abundantly with water upon disposal if reagents are spilled.
- R1/dilution liquid contains detergents that are irritating to eyes and skin. Flush with water if spilled or in contact with skin.

TEST CHARACTERISTICS

Analytical specificity

Monoclonal antibodies specific to human CRP are used in the test. No other human blood components are found to cross react with CRP in the NycoCard CRP System.

Standardization

NycoCard CRP is calibrated against ERM[®]-DA/IFCC reference preparation.

Measuring range

NycoCard READER II, SN ≥ 54170:

Whole blood samples: 8-200 mg/L

Serum or plasma samples: 5-120 mg/L

Measuring interval: 1 mg/L

Reference range

< 5 mg/L¹

Precision

A precision study was run according to CLSI guideline EP5-A2. Within-device (total) precision was determined for a plasma sample and Afinion CRP Control levels C I and C II

for 20 days, each sample assayed in duplicate, twice a day. An EDTA-whole blood sample was assayed for 4 days, each sample assayed in 20 replicates per day. The data are summarized in the table below.

Sample	N (days)	Mean CRP (mg/L)	Within-device precision (total) CV (%)
Serum	20	97	5.2
Control C I	20	22	6.4
Control C II	20	69	6.6
EDTA-whole blood	4	51	3.7

Limitations of the test

Heparin, citrate or EDTA, used as anticoagulant, does not interfere with the test. Increased levels of bilirubin, lipids or rheumatoid factor (RF) have none or negligible effect on the test results. Very high leukocytes levels (> 40x10⁹/L) in the sample might cause falsely elevated results (see "Trouble shooting"). Corrections should be made for hematocrit deviating from 40% (see "Effect of hematocrit").

STABILITY AND STORAGE

Unopened kits

The expiry dates apply to kit components stored in sealed, original containers at 2-8°C. Avoid direct sunlight or exposure to temperatures above 25°C. Do not freeze.

Opened kits

Avoid direct sunlight or exposure to temperatures above 25°C.

TD/Test device: Stable until expiry date at 2-8°C or 6 weeks at 15-25°C (refrigerate at the end of the working day). The sealing foil should be removed immediately before use.

R1/Dilution liquid: Stable until expiry date, refrigerated or at room temperature (2-25°C).

R2/Conjugate: Stable until expiry date at 2-8°C or 6 weeks at 15-25°C (refrigerate at the end of the working day).

R3/Washing solution: Stable until expiry date at 2-8°C or 6 weeks at 15-25°C (refrigerate at the end of the working day).

Samples

Whole blood (with anticoagulant): Stable for 3 days at 2-8°C.

Serum or plasma: Stable for 3 days at 2-8°C. For prolonged storage, keep samples frozen (-20°C).

Sample diluted with R1: Analyse diluted sample within 1 hour (diluted EDTA blood or plasma within 15 minutes).

TEST PROCEDURE

Important procedural notes!

- Bring test tube(s) with R1/dilution liquid to room temperature (15-25°C) before use.
- R2/conjugate, R3/washing solution and TD/test device can be used cold or equilibrated to room temperature.
- Mark the R1 tube and the test device with suitable patient or control identification.
- Always use a new pipette tip.
- Do not touch the test membrane with the pipette tip.
- The procedure steps should be performed successively without any interruptions.
- Caps should be screwed tightly on after use.

Sample material

Capillary blood, venous blood with or without anticoagulant (heparin, citrate or EDTA), serum and plasma (heparin, citrate or EDTA) can be used.

Quality control

The Afinion CRP Control should be used to confirm the efficacy of the reagents and correct performance of the test. Test this control according to the same procedure as for a patient sample.

Note! Use the “CRP Serum/Plasma” menu for reading the result with Nycocard Reader II. The measured value should be within the acceptance limits stated on the vial label.

Test procedure (illustrations on page 78)

1 Dilute sample

Fill a 5 µL capillary with patient sample or Afinion CRP Control, and drop the capillary into the tube with R1/dilution liquid. Close the tube and mix well for 10 seconds.

Note! Avoid air bubbles in the capillary and excess sample on the outside of the capillary.

2 Apply sample

Apply 50 µL diluted sample or diluted Afinion CRP Control to the TD/test device. Allow the sample to soak into the membrane (approx. 30 seconds).

Note! Avoid air bubbles on the membrane. Do not touch the membrane with the pipette tip.

3 Apply R2/conjugate

Apply one drop R2/conjugate to the TD/test device. Allow the reagent to soak into the membrane (approx. 30 seconds).

Note! The droplet bottle should be held vertically, about 1 cm above the membrane.

4 Apply R3/washing solution

Apply one drop R3/washing solution to the TD/test device. Allow the reagent to soak into the membrane (approx. 20 seconds).

Note! The droplet bottle should be held vertically, about 1 cm above the membrane.

5 Read the result

Read the result within 5 minutes using the Nycocard Reader II. Follow the Reader II User Instruction Manual.

Note! Use the “CRP Whole Blood” menu for reading whole blood samples, and the “CRP Serum/Plasma” menu for reading serum or plasma samples and the Afinion CRP Control.

INTERPRETATION OF RESULTS

Interpret the NycoCard CRP Test results with careful consideration to the patient's medical history, clinical examinations and other laboratory results. If the test result is questionable or if clinical signs and symptoms appear inconsistent with the test result, re-test the sample or confirm the result using another method. Analyse control materials frequently to verify the performance of the NycoCard Reader II Test System.

Effect of hematocrit

The "CRP Whole Blood" menu of the NycoCard Reader II is calibrated to read the CRP serum concentration of a blood sample with a hematocrit of 40%. If the actual hematocrit value deviates from 40%, the result should be corrected by multiplying with the respective factor in the table:

Hct%	Factor	Hct%	Factor
20-29	0.8	56-58	1.4
30-36	0.9	59-61	1.5
37-42	1.0	62-63	1.6
43-47	1.1	64-65	1.7
48-51	1.2	66-67	1.8
52-55	1.3	68-69	1.9

Reference range,

Hct: Women: 35-44%. Men: 39-48%.

TROUBLE SHOOTING

Possible cause of unexpected low result	Corrective action
Incorrect sample volume.	Retest the sample. Make sure that the capillary is filled completely. Avoid air bubbles.
Incorrect volume of diluted sample applied.	Check the pipette (change the Mini-Pet once a year). Make sure that the pipette tips fit well on the pipette. Retest the sample.
Too long delay (> 2 minutes) between steps 2 and 3.	Retest the sample. Steps 2 to 4 should be run without any interruptions.
Possible cause of unexpected high result	Corrective action
Very high leukocytes count (> 40x10 ⁹ /L) might give turbid sample dilutions. The time before the sample is absorbed into the test device will increase due to a partially blocked membrane.	Clarify the diluted sample by centrifugation. Retest the sample. Alternatively, use serum or plasma.
Incorrect sample volume.	Retest the sample. Avoid excess sample on the outside of the capillary.
Incorrect volume of diluted sample applied.	Check the pipette (change the Mini-Pet once a year). Make sure that the pipette tips fit well on the pipette. Retest the sample.
A reagent has been applied to the TD/test device before the preceding reagent has soaked completely into the membrane.	Retest the sample. Make sure that each reagent soaks completely into the membrane before applying the next reagent.
Too long delay (> 2 minutes) between steps 3 and 4.	Retest the sample. Steps 2 to 4 should be run without any interruptions.

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© 2019 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.