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## SARSTEDT Safety-Lanzette

<b>DE</b>	Gebrauchshinweis – SARSTEDT Safety-Lanzette	1 - 4
<b>EN</b>	Instructions for use – SARSTEDT safety lancet	5 - 7
<b>BG</b>	Инструкции за употреба – SARSTEDT обезопасена ланцета	8 - 10
<b>CS</b>	Návod k obsluze – Bezpečnostní lanceta SARSTEDT	11 - 13
<b>DA</b>	Brugsanvisning – SARSTEDT sikkerhedslancet	14 - 16
<b>EL</b>	Οδηγίες χρήσεως – Σκαρφιστήρας ασφαλείας SARSTEDT	17 - 19
<b>ES</b>	Instrucciones de uso – Lanceta Safety SARSTEDT	20 - 22
<b>FR</b>	Mode d'emploi – Lancette de sécurité SARSTEDT	23 - 25
<b>HR</b>	Uputa za upotrebu – SARSTEDT sigurnosna lanceta	26 - 28
<b>HU</b>	Használati utasítás – SARSTEDT biztonságos lándzsa	29 - 31
<b>IT</b>	Istruzioni d'uso – Lancetta Safety SARSTEDT	32 - 34
<b>KO</b>	사용 설명서 – SARSTEDT 안전 랜셋	35 - 37
<b>LV</b>	Lietošanas norādes – SARSTEDT drošības lancete	38 - 40
<b>NL</b>	Gebruiksaanwijzing – SARSTEDT Safety-lancet	41 - 43
<b>NO</b>	Bruksanvisning – SARSTEDT Safety-lansett	44 - 46
<b>PL</b>	Instrukcja obsługi – Bezpieczny lancet SARSTEDT	47 - 49
<b>PT</b>	Instruções de utilização – Lanceta de segurança SARSTEDT	50 - 52
<b>RO</b>	Instrucțiuni de utilizare – Lansetă de siguranță SARSTEDT	53 - 55
<b>RU</b>	Инструкция по применению – Безопасный ланцет SARSTEDT	56 - 58
<b>SK</b>	Návod na Použitie – SARSTEDT Bezpečnostná lanceta	59 - 61
<b>SL</b>	Navodila za uporabo – Varna lanceta SARSTEDT	62 - 64
<b>SV</b>	Bruksanvisning – SARSTEDT Safety-lansett	65 - 67
<b>TH</b>	คำแนะนำในการใช้ งาน – ใบมีดโรยตัวของ SARSTEDT	68 - 70
<b>TR</b>	Kullanım için talimatlar – SARSTEDT Güvenlikli lanset	71 - 73
<b>ZH</b>	使用说明 – SARSTEDT安全采血针	74 - 76

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### Normen zu Sterilisation:

DIN EN ISO 11137-1: Sterilization of health care products - Radiation - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.

DIN EN ISO 11137-2: Sterilization of health care products - Radiation - Part 2: Establishing the sterilization dose

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DIN EN ISO 11607-1: Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems

DIN EN ISO 11607-2: Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes

### Systembezogene Normen:

DIN EN ISO 13485: Medical devices - Quality management systems – Requirements for regulatory purposes

DIN EN ISO 14971: Medical devices - Application of risk management to medical devices

CLSI\* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.

CLSI\* GP42 "Collection of Capillary Blood Specimens", Approved Standard.

\*CLSI (Clinical and Laboratory Standards Institute)

MDSAP: *Medical Device Single Audit Program* gemäß ISO 13485

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




## Intended use

The safety lancet is used for the puncture to collect capillary blood. The safety lancet, with a prepared sample tube, e.g. micro sample tube or a POC test, is a system for collecting capillary blood. The products are intended for use in a professional environment by qualified medical personnel.

## Product description

The safety lancet is a sterile product for single use. The safety lancet consists of a plastic casing, a disposable needle/blade of stainless steel, a release button in various colours, steel springs and a safety plunger.

## Safety lancet product line

Image	Version	Needle size	Puncture depth	Blood volume approx.*	
				After 20 sec.	After 120 sec.
	Mini	28 G	1.6 mm	50 µl	150 µl
	Normal	21 G	1.8 mm	100 µl	500 µl
	Extra	18 G	1.8 mm	150 µl	750 µl
	Super	Blade 1.5 mm	1.6 mm	150 µl	900 µl
	Neonatal	Blade 1.5 mm	1.2 mm	100 µl	450 µl

\*According to a study with 30 adult subjects.

We would like to point out that the selection of the correct version depends on various factors such as skin type, skin blood flow, puncture site and the correct placement of the safety lancet on the skin. If in doubt, testing several different safety lancets is recommended before routine use. The information presented here is merely a recommendation and is not a substitute for medical, scientific or technical advice.

## Safety instructions and warnings

1. General precautions: Wear gloves and other general personal protective equipment to protect yourself from blood and possible exposure to pathogens transmissible by blood.
2. Treat all biological specimens and sharp blood collection utensils (needles, adapters) according to the guidelines and procedures in your institution. In case of direct contact to biological specimens or a needle prick injury, consult a doctor, as HIV, HCV, or HBV or other infectious diseases can be transmitted. The safety guidelines and procedures of your institution must be followed.
3. Dispose of the safety lancet in suitable containers for sharp objects.
4. Before using, check the safety plunger for damage and do not use if it is damaged.
5. Reusing the product can lead to infections, injuries or death.
6. Do not use the product after the expiry date. The expiry date of the product is the last day of the month and year indicated.

## Storage

Store the product at room temperature.

## Handling

**READ THIS DOCUMENT COMPLETELY BEFORE STARTING THE PUNCTURE.**

### Materials needed for collecting a specimen:

1. All necessary sample containers, e.g. micro sample tubes, labelled by size and additive as per request.
2. Gloves, lab coat, eye protection or other suitable protective gear to protect against pathogens transmitted by blood or potentially infectious material.
3. Labels for sample identification.
4. Safety lancet.
5. Disinfectant for cleaning the puncture site (follow the facility's guidelines for preparing the puncture site for taking blood samples). Do not use any cleaning materials containing alcohol if the specimens are to be used for a blood alcohol test.
6. Dry, germ-free swabs.
7. Plaster or dressing.
8. Sharps disposal container for the safe disposal of used material.




**NOTE: Follow the regulations of your institution for the order of collection.**

## Puncture with the safety lancet

### General instructions:

WEAR GLOVES TO MINIMISE THE RISK OF EXPOSURE DURING CAPILLARY BLOOD COLLECTION.

1. Select a suitable Microvette® for the required specimen.
2. Select a site for the puncture with a safety lancet. Please follow the guidelines of the institution for preparing the puncture site.

- 1 Twist the protective cap off.
- 2 Hold the safety lancet against the selected and disinfected puncture site. Press the release button.
- 3 Dispose of the safety lancet in a suitable container for sharp objects.

The different versions of the safety lancet are suitable as needed for collecting small or medium up to large volumes of blood (see: Product description "Safety lancet product range")

### Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Dispose of the safety lancet in a sharps container.
3. Disposable gloves prevent the risk of infection.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

## US standards / ISO standards in the latest version

### Standards for labelling:

DIN EN 1041: Information supplied by the manufacturer of medical devices

DIN EN ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements

### Standards for needles/cannulas:

DIN EN 10088-1, Table 2 - Stainless steels - Part 1: List of stainless steels

### Standards for sterilisation:

DIN EN ISO 11137-1: Sterilization of health care products - Radiation - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.

DIN EN ISO 11137-2: Sterilization of health care products - Radiation - Part 2: Establishing the sterilization dose

### Standards for packaging:

DIN EN ISO 11607-1: Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems

DIN EN ISO 11607-2: Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes

### System-related standards:

DIN EN ISO 13485: Medical devices - Quality management systems – Requirements for regulatory purposes

DIN EN ISO 14971: Medical devices - Application of risk management to medical devices

CLSI\* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.

CLSI\* GP42 "Collection of Capillary Blood Specimens", Approved Standard.

\*CLSI (Clinical and Laboratory Standards Institute)

MDSAP: *Medical Device Single Audit Program* in accordance with ISO 13485

MDR 93/42/EEC for medical devices

## Key for icons and labels:

Article number

Batch number

CE marking

Manufacturer

Sterilisation by irradiation

Risk of contamination if reused

Do not resterilise

Country of manufacturing

Use before

Follow the instructions for use

Store away from sunlight

Store in a dry place

Do not use if package is damaged

Medical device

Single sterile barrier system with external protective packaging

Technical modifications reserved

All serious incidents relating to the product shall be notified to the manufacturer and the competent national authority.

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EN 1041: Information supplied by the manufacturer of medical devices

EN ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements

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### Normy pro sterilizaci:

EN ISO 11137-1: Sterilization of health care products - Radiation - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.

EN ISO 11137-2: Sterilization of health care products - Radiation - Part 2: Establishing the sterilization dose

### Normy pro obal:

EN ISO 11607-1: Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems

EN ISO 11607-2: Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes

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EN ISO 13485: Medical devices - Quality management systems - Requirements for regulatory purposes

EN ISO 14971: Medical devices - Application of risk management to medical devices

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MDSAP: *Medical Device Single Audit Program* podle ISO 13485

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DIN EN 1041 : Information supplied by the manufacturer of medical devices

DIN EN ISO 15223-1 : Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements

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DIN EN ISO 11137-1 : Sterilization of health care products - Radiation - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.

DIN EN ISO 11137-2 : Sterilization of health care products - Radiation - Part 2: Establishing the sterilization dose

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DIN EN ISO 11607-1 : Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems

DIN EN ISO 11607-2 : Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes

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DIN EN ISO 13485 : Medical devices - Quality management systems – Requirements for regulatory purposes

DIN EN ISO 14971 : Medical devices - Application of risk management to medical devices

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DIN EN 1041: Information supplied by the manufacturer of medical devices  
 DIN EN ISO 15223-1: Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements

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DIN EN ISO 11137-1: Sterilization of health care products – Radiation – Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.  
 DIN EN ISO 11137-2: Sterilization of health care products – Radiation – Part 2: Establishing the sterilization dose

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DIN EN ISO 11607-1: Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems  
 DIN EN ISO 11607-2: Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes

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DIN EN ISO 13485: Medical devices – Quality management systems – Requirements for regulatory purposes  
 DIN EN ISO 14971: Medical devices – Application of risk management to medical devices  
 CLSI\* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.  
 CLSI\* GP42 "Collection of Capillary Blood Specimens", Approved Standard.  
 \*CLSI (Clinical and Laboratory Standards Institute)  
 MDSAP: ISO 13485에 따른 Medical Device Single Audit Program  
 의료 기기용 MDD 지침 93/42/EWG

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DIN EN 1041: Information supplied by the manufacturer of medical devices  
 DIN EN ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements

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DIN EN ISO 11137-1: Sterilization of health care products - Radiation - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.  
 DIN EN ISO 11137-2: Sterilization of health care products - Radiation - Part 2: Establishing the sterilization dose

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DIN EN ISO 11607-1: Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems  
 DIN EN ISO 11607-2: Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes

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DIN EN ISO 13485: Medical devices - Quality management systems – Requirements for regulatory purposes  
 DIN EN ISO 14971: Medical devices - Application of risk management to medical devices  
 CLSI\* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.  
 CLSI\* GP42 "Collection of Capillary Blood Specimens", Approved Standard.  
 \*CLSI (Clinical and Laboratory Standards Institute)  
 MDSAP: Medical Device Single Audit Program ตามมาตรฐาน ISO 13485  
 ระเบียบข้อบังคับ MDD 93/42/EWG สำหรับผลิตภัณฑ์ทางการแพทย์

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DIN EN 1041:Information supplied by the manufacturer of medical devices

DIN EN ISO 15223-1:Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1:General requirements

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DIN EN ISO 13485: Medical devices - Quality management systems - Requirements for regulatory purposes

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CLSI\* GP39 " Tubes and Additives for Venous and Capillary Blood Specimen Collection ", Approved Standard.

CLSI\* GP42 " Collection of Capillary Blood Specimens ", Approved Standard.

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MDSAP: 根据 Medical Device Single Audit Program 规定, 符合 ISO 13485

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