PERFORMANCES ANALYTIQUES DU TEST

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*NMT : niveau maximum testé

RÉFÉRENCES

- Cancer Research Campaign Factsheet No. 1, 1994.
- Cancer Facts and Figures**, American Cancer Society, 1996.
- Thrasher J., Crawford E.: Current Management of Invasive Metastatic Transitional Cell Carcinoma of the Bladder. **The Journal of Urology** 149:957-972, 1993.
- Murphy W.: Current Status of Urinary Cytology in the Evaluation of Bladder Neoplasms. **Human Pathology** 21:886-896, 1990.
- Umiker W.: Accuracy of Cytologic Diagnosis of Cancer of the Urinary Tract. **Symposium on Diagnostic Accuracy of Cytologic Technics** 8:186-193,1964.
- Badalament R.A., Hermansen D.K., Kimmel M., Gay H., Herr H.W., Fair W.R., Whitmore W.F., Jr., Melamed M.R.: The Sensitivity of Bladder Wash Flow Cytometry, Bladder Wash Cytology, and Voided Cytology in the Detection of Bladder Carcinoma. **Cancer** 60:1423-1427, 1987.
- Sarosdy M.F., Hudson M.A., et al: Improved Detection of Recurrent Bladder Cancer Using the Bard BTA *stat* Test. **Urology** 50, 349-353, 1997.
- Raitanen M.P., Marttila T., et al: The Bard BTA *stat* Test in Monitoring Bladder Cancer. **Journal of Urology** 159 (No. 5): 244, 1998.
- Leyh H., Treiber U.: European Multicenter Trial Compares the BTA *stat* Test to Urine Cytology in Patients Suspected of Having Bladder Cancer. **European Urology** 33 (Supplement 1): 140, 1998.
- Corey M., Kinders R., et al: Factor H Related Proteins Are Upregulated In Bladder Cancer. **Proceedings of the American Association for Cancer Research** 39: 263, 1998.
- Kinders R., Jones T., et al: Complement Factor H or a Related Protein Is a Marker for Transitional Cell Cancer of the Bladder. **Clinical Cancer Research** 4: 2511-2520, 1998.
- Austyn J. M., Wood K. J.: **Principles of Cellular and Molecular Immunology**. Oxford University Press p. 522-554, 1993.

Distributed by:

 Polymedco, Inc. European Authorized Representative

Obelis s.a. (O.E.A.R.C.)
Bd Général Wahnis 53
1030 Brussels
Belgium
Tel: +(32) 2 732-59-54
Fax: +(32) 2 732-60-03

Manufactured for:
Polymedco, Inc.
510 Furnace Dock Road
Cortlandt Manor, NY 10567
U.S.A.

FM 77123

© 2012
BTA *stat* is a registered trademark of Polymedco, Inc.
U.S. Patent No. 2,111,161 and 6,017703
Enregistrement AFSSAPS n° R51282

Issued 12/12
PN 52401-00

BTA *stat* TEST

Test Kit No. HI-308-37050601 (10 Tests)

Test Kit No. HI-308-37050602 (5 Tests)

A Rapid Test For the Qualitative Detection of
Bladder Tumor Antigen in Human Urine

Made in the United States.

Multi-language labeling for export only.

For professional use only.

INTENDED USE

The BTA *stat* test is an *in vitro* diagnostic assay indicated for the qualitative detection of bladder tumor antigen in human urine. This test is intended for use as an aid in the diagnosis and management of bladder cancer patients.

SUMMARY AND EXPLANATION OF THE TEST

Bladder Cancer

Bladder cancer is the fifth most common form of cancer in men and the eleventh most common in women in Europe.¹ The overall incidence is approximately 19 per 100 000. Similar to the prevalence in Europe, bladder cancer ranks fourth among men and ninth among women in the United States.² Approximately 75 to 85% of these patients present with transitional cell carcinoma (TCC) confined to the superficial mucosa of the bladder.³ The risk of recurrence in these patients is 75%. Patients with previous diagnosis of bladder cancer have been routinely followed for recurrence by urine cytology and cystoscopy. Both methods have their limitations.

Cystoscopy is considered the diagnostic standard for sensitivity and specificity when a biopsy is not obtained. This method is an invasive procedure associated with patient discomfort, expense and is limited to diagnosis of those tumors that can be visualized.⁴

Voided urine cytology (VUC), or the examination of urinary sediment for cancer cells has several factors that contribute to suboptimal results. Urothelial cells require about a year to replicate, so few are available for examination in any particular sample.⁵ Exfoliated cells enter a hostile environment of high acidity and low osmolality which may obscure essential diagnostic features. Standards for specimen collection, preservation, processing and interpretation have not been widely accepted. Routine cytology, as a monitoring tool exhibits variable sensitivity depending on the tumor stage and grade with lowest sensitivity reported for early stage disease.^{5,6} In addition, the best quality results are obtained from examination of samples collected under specific void procedures.⁵

The diagnosis and management of bladder cancer patients could be improved with a rapid, simple, urine test that could be performed at point of care or in the laboratory. Recent studies have shown that the BTA *stat* test can be extremely useful in this regard.^{7,9} The BTA *stat* test, which qualitatively detects bladder tumor antigen in urine, is a single-step, antibody-based test which can be performed in only 5 minutes with no pretreatment of the urine sample.

Bladder Tumor Antigen

The monoclonal antibodies used in the BTA *stat* test were generated against urine components from bladder cancer patients who were positive by histology. The bladder tumor antigen was identified as human complement factor H-related protein (hCFHrp) similar in composition, structure, and function to human complement factor H (hCFH).^{10,11} hCFH is also recognized by the monoclonal antibodies used in the BTA *stat* test. In cell culture, hCFHrp was shown to be produced by several human bladder cancer cell lines but not by normal epithelial cell lines.^{10,11}

hCFH plays a key inhibitory role in the control of the alternative complement pathway that functions to lyse cells recognized as foreign to the host. By interaction with complement factor C3b, hCFH serves to inhibit the formation of a membrane attack complex, thereby preventing cell lysis.¹² *In vitro*, bladder tumor antigen interrupts the complement cascade and protects cells from lysis by complement. Production of bladder tumor antigen may confer a selective growth advantage to cancer cells *in vivo* by allowing the cells to evade the host immune system.

PRINCIPLE OF THE PROCEDURE

The BTA *stat* rapid test for bladder tumor antigen is an immunochromatographic assay utilizing monoclonal antibodies to specifically detect the presence of bladder tumor antigen in urine. Patient urine is added to the sample well and allowed to react with a colloidal gold-conjugated antibody. If the antigen is present in the sample, an antigen-conjugate complex is formed. The reaction mixture flows through the membrane which contains zones of immobilized antibodies. In the Patient (P) zone, antigen-conjugate complexes are captured by a second antibody, forming a visible line. In the absence of the antigen in the patient urine, no visible line will form. The Control (✓) zone contains an immobilized reagent which will capture the conjugate independently of the presence or absence of the antigen, thereby always producing a line. This procedural control assures the operator that each device is working properly.

CONTRAINDICATIONS

- Do not use beyond the labeled expiration date.
- Do not reuse disposable test devices. Discard after single use.
- Do not use if pouch is damaged or opened.
- Do not touch the membrane located within the windows.

WARNINGS

- For *in vitro* diagnostic use.

PRECAUTIONS

- To avoid cross-contamination of samples, use a new dropper (provided with each device) for each patient urine.
- Treat urine samples and used devices as if they are potentially infectious.

STORAGE AND STABILITY

- Store the BTA *stat* Test Kit at 2 - 30°C. The Test Kit is stable when stored at these temperatures until the expiration date printed on the box label.
- Do not freeze the BTA *stat* Test Kit.

SPECIMEN COLLECTION, STORAGE AND PREPARATION

Voided urine or urine from a catheterized patient is required for the BTA *stat* test. Bladder barbotage specimens, serum, plasma or whole blood should not be used. Urine should be collected without preservatives or fixatives in a clean urine cup and labeled appropriately. If urine is to be used for other tests, remove an aliquot of specimen (a minimum of 2 mL) for this test to avoid contamination. Mix urine before testing. Urine samples may be stored at room temperature for up to 48 hours after collection. If the urine sample is not tested within 48 hours, it should be refrigerated at 2 - 8°C for up to 7 days. If the refrigerated urine sample is not tested within 7 days, it should be stored at or below -20°C until tested. Frozen samples may be stored indefinitely. Avoid repeated freezing and thawing.

- Do not use paper or foam cups for urine specimen collection or storage.
- The effect of radiation therapy or systemic chemotherapy on the BTA *stat* test is unknown.
- Patients receiving intravesical agents, such as BCG, mitomycin C, Thiotepa, bropiramine (investigational) or interferon (investigational), should be tested at the discretion of the physician.
- The effects of experimental drugs on the BTA *stat* test are unknown.
- Some patients with benign renal disease (such as stones and nephritis) or renal cancer may test positive in the BTA *stat* test.
- For trauma to the bladder or urinary tract due to surgery, biopsy, etc. the physician should allow ample time for trauma recovery before using the test.

CONTENTS OF KIT

BTA *stat* Test Kit (Cat. No. HI-308-37050601 / Cat. No. HI-308-37050602)

Components:


- 10 Foil Packages (HI-308-37050601) / 5 Foil Packages (HI-308-37050602)
Each package contains -
 - 1 BTA *stat* device. Each device incorporates a membrane with capture antibody and second conjugated antibody in a protein matrix containing sodium azide.
 - 1 Disposable dropper
 - 1 Disposable desiccant pouch
- 1 Package Insert.

MATERIALS REQUIRED BUT NOT PROVIDED

- Timer
- Urine collection container

PATIENT TEST PROCEDURE


1. Bring test materials and patient urine sample to room temperature. Gently swirl patient's urine sample.
2. Remove the test device and dropper from foil package. Throw away small desiccant pouch. Place the device on a clean, well-lit, flat surface and label with the patient's identification.
3. Set timer. Fill the dropper provided with the patient's urine sample and hold it upright above the sample well as shown.
4. Allow 3 **FULL** drops (without air bubbles) to fall into the sample well. Start timer. Throw dropper away.
5. When timer reaches 5 minutes, **read results within 1 minute**. Read results as shown under "Interpretation of Results."


Read at 5 minutes but **NO LATER THAN 6 MINUTES**. Test result is not valid if read after 6 minutes.

INTERPRETATION OF RESULTS


Positive Result: Carefully look at Patient (P) window. **ANY** pink or reddish-brown colored line, **NO MATTER HOW FAINT**, in the Patient (P) window is a positive result.

Negative Result: Carefully look at Patient (P) window. No colored line in the Patient (P) window is a negative result.

Invalid Test Result: If no line appears in the Control (✓) window, the test is invalid and must be repeated with a new device. The most common reason for an invalid test result is failure to add exactly 3 FULL drops.

QUALITY CONTROL

Good laboratory practices recommend the use of appropriate controls.

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LIMITATIONS

Results of the BTA stat test should not be interpreted as absolute evidence for the presence or absence of TCC of the bladder. Elevated levels of bladder tumor antigen have been observed in patients with renal stones, nephritis, renal cancer and recent trauma to the bladder or urinary tract.

The result from the BTA stat test should be used only in conjunction with information available from the clinical evaluation of the patient and other diagnostic procedures.

CLINICAL PERFORMANCE

CLINICAL SPECIFICITY

BTA stat test specificity was studied using urine samples from 555 non-bladder cancer patients. Samples were collected from diverse geographic locations and stored frozen until tested. The average age was 55 years and 48% were males. The normal healthy population consisted of smokers and non-smokers. The non-genitourinary diseases included diabetes, arthritis, lupus and other collagen degenerative diseases, and malignancies included mammary, pulmonary, gastrointestinal carcinomas, leukemia and lymphomas. The non-bladder genitourinary (GU) cancers category consisted of prostate, renal cell, endometrial, ovarian and other GU carcinomas. The GU disease category consisted of patients with prostatitis, urethritis, renal stones and disease, urinary tract infections, incontinence, sexually transmitted diseases and other disorders.

The results indicated that in healthy individuals and individuals without GU diseases and malignancies, the BTA stat test negative rate was 95% and 93%, respectively. Positive BTA stat test results may occur in patients with renal disease such as stones and nephritis and patients with renal cancer (including upper tract TCC).

The BTA stat test specificity in 107 patients with a history of bladder cancer and no evidence of disease was 70%.

Expected results may vary depending on the patient population tested.

Table I. BTA stat TEST CLINICAL SPECIFICITY RESULTS

PATIENT TYPE	NUMBER OF SUBJECTS	TEST NEGATIVE (%)
Healthy Subjects	167	95
Non-Genitourinary Diseases and Cancers	105	93
Non-Genitourinary Diseases	52	98
Non-Genitourinary Cancers	53	89
Genitourinary Diseases	152	72
Benign Prostatic Hyperplasia	26	88
Benign Renal Disease	32	50
Misc. Genitourinary Disease	94	76
Genitourinary Cancers	77	73
Prostate Cancers	45	78
Renal Cancers	7	29
Other Cancers	25	76
Genitourinary Trauma	54	33

CLINICAL SENSITIVITY

BTA stat test sensitivity results were determined using urine samples from 220 patients with biopsy proven bladder tumors. Samples were collected from diverse geographic locations and stored frozen until tested. The average age was 68 years and 79% were males. Results are presented below by stage and by grade of the tumor.

Table II. BTA stat TEST CLINICAL SENSITIVITY BY STAGE AND GRADE*

STAGE	N	SENSITIVITY (%)
Ta	111	51
T1	38	90
≥T2	50	88
Tis	18	61
GRADE	N	SENSITIVITY (%)
1	57	42
2	56	66
3	95	83

*3 patients without stage and 12 without grade determinations.

A subset of the patients studied (131) also had voided urine cytology (VUC) performed on the same sample as the BTA stat test. The BTA stat test was shown to be more sensitive than VUC in all categories except for Tis.

Table III. BTA stat TEST AND VUC CLINICAL SENSITIVITIES

STAGE	N	SENSITIVITY BTA stat (%)	SENSITIVITY VUC (%)	SENSITIVITY BTA stat + VUC (%)
Ta	73	45	7	49
T1	27	85	41	85
≥T2	16	75	38	81
Tis	15	53	60	80

ANALYTICAL PERFORMANCE

HIGH DOSE HOOK EFFECT

High dose hook (prozone) effect tests were conducted to determine if the BTA *stat* test is free from interference from high concentration positive patient samples. Results showed that there was no prozone effect up to 12400 U/mL bladder tumor antigen in a patient's urine sample, which was the highest concentration available for testing.

REPRODUCIBILITY

Inter-assay reproducibility studies were conducted using two lots of BTA *stat* devices and lot-to-lot reproducibility studies were conducted using three lots of BTA *stat* devices to test 10 replicates of 4 blinded samples per day for 5 days using three independent readers for each lot of devices. The results of the studies are shown in Table IV below.

Table IV. REPRODUCIBILITY OF THE BTA *stat* TEST

INTER-ASSAY			Lot: Dev 3		Lot: Dev 4	
Sample	Expected Result	N	No. Pos	No. Neg	No. Pos	No. Neg
A	N	50*	0	50	0	50
C	P	50	50	0	50	0
E	P	50	50	0	50	0
F	P	50	50	0	50	0

LOT-TO-LOT			Lot: Dev 3		Lot: Dev 4		Lot: Pilot 1	
Sample	Expected Result	N	No. Pos	No. Neg	No. Pos	No. Neg	No. Pos	No. Neg
A	N	50*	0	50	0	50	0	50
C	P	50	50	0	50	0	49	1
E	P	50	50	0	50	0	49	1
F	P	50	50	0	50	0	50	0

*10 per day, 5 days

Sample C contains an amount of bladder tumor antigen near the limit of detection of the BTA *stat* test.

LIMIT OF DETECTION

The limit of detection was established at 22,4 U/mL.

INTERFERING SUBSTANCES

Normal and TCC positive urine pools containing the substances listed below were tested in the BTA *stat* test. These substances did not show interference in the BTA *stat* test at the levels indicated.

SUBSTANCE	HIGHEST LEVEL TESTED WITH NO INTERFERENCE	LEVEL AT WHICH SUBSTANCE INTERFERED
Possible Urine Constituents		
Hemoglobin	100 mg/dL	No interference at MLT*
White Blood Cells	10 ⁶ cells/mL	No interference at MLT
Red Blood Cells	10 ⁶ cells/mL	No interference at MLT
Albumin	1 g/dL	No interference at MLT
Bilirubin (unconjugated)	0,4 mg/dL	0,8 mg/dL
IgG	10 mg/dL	No interference at MLT
Uric Acid	250 mg/dL	No interference at MLT
Ascorbic Acid	5 g/dL	No interference at MLT
Therapeutic Agents		
Ampicillin	600 mg/dL	No interference at MLT
Doxorubicin-HCl	10 mg/dL	No interference at MLT
Mitomycin C	10 mg/dL	No interference at MLT
Nitrofurantoin	50 mg/dL	No interference at MLT
Phenazopyridine-HCl	80 mg/dL	100 mg/dL
Thiotepa	10 mg/dL	No interference at MLT
Trimethoprim	50 mg/dL	No interference at MLT
Bacillus Calmette Guerin	20 mg/dL	No interference at MLT
Finasteride	2,5 mg/dL	No interference at MLT
Flutamide	100 mg/dL	No interference at MLT
Ioversol, 74% (imaging contrast agent)	1%	5%
Urised	17,5 mg/dL	35 mg/dL

*MLT - Maximum Level Tested

REFERENCES

- Cancer Research Campaign Factsheet No. 1, 1994.
- Cancer Facts and Figures**, American Cancer Society, 1996.
- Thrasher J., Crawford E.: Current Management of Invasive Metastatic Transitional Cell Carcinoma of the Bladder. **The Journal of Urology** 149:957-972, 1993.
- Murphy W.: Current Status of Urinary Cytology in the Evaluation of Bladder Neoplasms. **Human Pathology** 21:886-896, 1990.
- Umiker W.: Accuracy of Cytologic Diagnosis of Cancer of the Urinary Tract. **Symposium on Diagnostic Accuracy of Cytologic Technics** 8:186-193,1964.
- Badalament R.A., Hermansen D.K., Kimmel M., Gay H., Herr H.W., Fair W.R., Whitmore W.F., Jr., Melamed M.R.: The Sensitivity of Bladder Wash Flow Cytometry, Bladder Wash Cytology, and Voided Cytology in the Detection of Bladder Carcinoma. **Cancer** 60:1423 - 1427, 1987.
- Sarosdy M.F., Hudson M.A., et al: Improved Detection of Recurrent Bladder Cancer Using the Bard BTA *stat* Test. **Urology** 50, 349 - 353, 1997.
- Raitanen M.P., Marttila T., et al: The Bard BTA *stat* Test in Monitoring Bladder Cancer. **Journal of Urology** 159 (No. 5): 244, 1998.
- Leyh H., Treiber U.: European Multicenter Trial Compares the BTA *stat* Test to Urine Cytology in Patients Suspected of Having Bladder Cancer. **European Urology** 33 (Supplement 1): 140, 1998.
- Corey M., Kinders R., et al: Factor H Related Proteins Are Upregulated In Bladder Cancer. **Proceedings of the American Association for Cancer Research** 39: 263, 1998.
- Kinders R., Jones T., et al: Complement Factor H or a Related Protein Is a Marker for Transitional Cell Cancer of the Bladder. **Clinical Cancer Research** 4: 2511-2520, 1998.
- Austyn J. M., Wood K. J.: **Principles of Cellular and Molecular Immunology**. Oxford University Press p. 522-554, 1993.

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LITERATURHINWEIS

- Cancer Research Campaign Factsheet No. 1, 1994.
- Cancer Facts and Figures**, American Cancer Society, 1996.
- Thrasher J., Crawford E.: Current Management of Invasive Metastatic Transitional Cell Carcinoma of the Bladder. **The Journal of Urology** 149:957-972, 1993.
- Murphy W.: Current Status of Urinary Cytology in the Evaluation of Bladder Neoplasms. **Human Pathology** 21:886-896, 1990.
- Umiker W.: Accuracy of Cytologic Diagnosis of Cancer of the Urinary Tract. **Symposium on Diagnostic Accuracy of Cytologic Technics** 8:186-193,1964.
- Badalament R.A., Hermansen D.K., Kimmel M., Gay H., Herr H.W., Fair W.R., Whitmore W.F., Jr., Melamed M.R.: The Sensitivity of Bladder Wash Flow Cytometry, Bladder Wash Cytology, and Voided Cytology in the Detection of Bladder Carcinoma. **Cancer** 60:1423- 1427, 1987.
- Sarosdy M.F., Hudson M.A., et al: Improved Detection of Recurrent Bladder Cancer Using the Bard BTA *stat* Test. **Urology** 50, 349- 353, 1997.
- Raitanen M.P., Marttila T., et al: The Bard BTA *stat* Test in Monitoring Bladder Cancer. **Journal of Urology** 159 (No. 5): 244, 1998.
- Leyh H., Treiber U.: European Multicenter Trial Compares the BTA *stat* Test to Urine Cytology in Patients Suspected of Having Bladder Cancer. **European Urology** 33 (Supplement 1): 140, 1998.
- Corey M., Kinders R., et al: Factor H Related Proteins Are Upregulated In Bladder Cancer. **Proceedings of the American Association for Cancer Research** 39: 263, 1998.
- Kinders R., Jones T., et al: Complement Factor H or a Related Protein Is a Marker for Transitional Cell Cancer of the Bladder. **Clinical Cancer Research** 4: 2511-2520, 1998.
- Austyn J. M., Wood K. J.: **Principles of Cellular and Molecular Immunology**. Oxford University Press p. 522-554, 1993.

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REFERENCIAS

- Cancer Research Campaign Factsheet No. 1, 1994.
- Cancer Facts and Figures**, American Cancer Society, 1996.
- Thrasher J., Crawford E.: Current Management of Invasive Metastatic Transitional Cell Carcinoma of the Bladder. **The Journal of Urology** 149:957-972, 1993.
- Murphy W.: Current Status of Urinary Cytology in the Evaluation of Bladder Neoplasms. **Human Pathology** 21:886-896, 1990.
- Umiker W.: Accuracy of Cytologic Diagnosis of Cancer of the Urinary Tract. **Symposium on Diagnostic Accuracy of Cytologic Technics** 8:186-193,1964.
- Badalament R.A., Hermansen D.K., Kimmel M., Gay H., Herr H.W., Fair W.R., Whitmore W.F.,Jr., Melamed M.R.: The Sensitivity of Bladder Wash Flow Cytometry, Bladder Wash Cytology, and Voided Cytology in the Detection of Bladder Carcinoma. **Cancer** 60:1423- 1427, 1987.
- Sarosdy M.F., Hudson M.A., et al: Improved Detection of Recurrent Bladder Cancer Using the Bard *BTA stat* Test. **Urology** 50, 349- 353, 1997.
- Raitanen M.P., Marttila T., et al: The Bard *BTA stat* Test in Monitoring Bladder Cancer. **Journal of Urology** 159 (No. 5): 244, 1998.
- Leyh H., Treiber U.: European Multicenter Trial Compares the *BTA stat* Test to Urine Cytology in Patients Suspected of Having Bladder Cancer. **Europeen Urology** 33 (Supplement 1): 140, 1998.
- Corey M., Kinders R., et al: Factor H Related Proteins Are Upregulated In Bladder Cancer. **Proceedings of the American Association for Cancer Research** 39: 263, 1998.
- Kinders R., Jones T., et al: Complement Factor H or a Related Protein Is a Marker for Transitional Cell Cancer of the Bladder. **Clinical Cancer Research** 4: 2511-2520, 1998.
- Austyn J. M., Wood K. J.: **Principles of Cellular and Molecular Immunology**. Oxford University Press p. 522-554, 1993.

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